



JWE GmbH

## Staining test apparatus



### Function

This device was developed to coat human or animal teeth with stain. The teeth are embedded in a resin. The samples are immersed alternately in a liquid of coffee, tea and gastric juice and dried. For drying, we offer the opportunity with heating lamps or a heating chamber. In addition there is the unit with a heater, for heating the liquid. The unit can hold up to 100 products. The samples will be cleaned with the brushing machine (see Brushing machine for the removal of dental plaque). The plaque or stain removal can be measured with an image analyzing system (see image capture of the remote plaque / Image analysis of Stain Removal).

The sample holder and the tray can be easily removed for cleaning and assembly by opening clamp levers.

We are happy to customize the test equipment to the geometric conditions of your required tensile testing machine.

Please contact us.



### Specification

Dimension (W,D,H) ca.	400, 300, 500 mm
Weight	20 kg
Power	2x250 W
Revolutions	1-5 rpm
Cycles	1 bis 9999
Teststations	4 x 25
Power Supply	115/230 (50/60 Hz)

### Order Information

Device with Heating Lamp	312-010
Device with Heating Chamber	312-020
Basic unit with liquid heating	312-030

### Scope of Delivery

1	Device
1	Operating Instruction
1	Power Cord CEE 16A (5 m)

### Conditions of Delivery

Delivery:	EXW
Shipping time:	by arrangement

### Contakt

JWE GmbH  
Tel. 0049 (0)6173 3948596  
Fax 032 22 1455 209  
Email [info@jw-e.net](mailto:info@jw-e.net)  
Web [www.jw-e.net](http://www.jw-e.net)

Technical changes reserved.  
The illustrations may differ from the actual equipment.  
All data without obligation and not-binding.